

## Spritzflaschen, IDL Wash Bottles

**Enghals-Spritzflasche LDPE**, rund, naturfarben, auslaufsicher mit beigefügtem Spritzverschluss LDPE, Steigrohr naturfarben.

| Volumen<br>ml | Verschluss-Ø<br>mm | Gewinde<br>A.-Ø<br>mm | Höhe<br>mm | Best.-Nr.          |
|---------------|--------------------|-----------------------|------------|--------------------|
| 100           | 18                 | 18                    | 105        | <b>10 00 80100</b> |
| 200           | 18                 | 18                    | 133        | <b>10 00 80200</b> |
| 250           | 25                 | 25                    | 140        | <b>10 00 80250</b> |
| 500           | 25                 | 25                    | 177        | <b>10 00 80500</b> |
| 1000          | 28                 | 28                    | 223,5      | <b>10 00 81000</b> |
| 2000          | 28                 | 28                    | 264        | <b>10 00 82000</b> |

